

# Update 3.2.0 qubeCLOUD



## Release Notes Version 3.2.0

### Allgemein

- Optimierter Login Dialog
- Optimierter Lizenzaktivierungsprozess und Darstellung der verfügbaren Simulationsmonitore
- Neues Mediacenter verfügbar

The screenshot displays the SkillQube SimConfig interface. The top navigation bar includes links for HOME, MY ORGANISATION, SIMCONFIG, EDUCATION CENTER, SCENARIOS, CUSTOMERS, QUBEPLAN, and ADMINISTRATION. The SIMCONFIG section is active, and the MEDIA tab is selected. Below the tabs, there is a 'Media' section with a 'CREATE' button and a search bar. A table lists the media files with columns for Name, User, Date, Scenario count, and Actions. The table contains 10 entries, including AHA ACLS, 04\_HWS Fraktur\_Horizont 1.2.png, 05\_HWS Fraktur\_Horizont 2.2.png, 09z.pdf, AHA 2020 - Brachycuria.pdf, AHA 2020 - Opliccos.pdf, AHA 2020 - Plano cardiaco.pdf, AHA 2020 - Resuscitation.pdf, AHA 2020 - Taquicardia.pdf, and AHA Algorithmus LVAD.pdf. The bottom of the page shows pagination controls: Items per page: 10, 1-10 of 60, and navigation arrows.

Name	User	Date	Scenario count	Actions
AHA ACLS	-	12/16/2025	-	[edit] [delete]
04_HWS Fraktur_Horizont 1.2.png	-	9/11/2025	2	[edit] [delete] [download]
05_HWS Fraktur_Horizont 2.2.png	-	9/11/2025	2	[edit] [delete] [download]
09z.pdf	-	9/11/2025	4	[edit] [delete] [download]
AHA 2020 - Brachycuria.pdf	-	9/11/2025	32	[edit] [delete] [download]
AHA 2020 - Opliccos.pdf	-	9/11/2025	8	[edit] [delete] [download]
AHA 2020 - Plano cardiaco.pdf	-	9/11/2025	98	[edit] [delete] [download]
AHA 2020 - Resuscitation.pdf	-	9/11/2025	38	[edit] [delete] [download]
AHA 2020 - Taquicardia.pdf	-	9/11/2025	28	[edit] [delete] [download]
AHA Algorithmus LVAD.pdf	-	9/11/2025	6	[edit] [delete] [download]



## SimKonfig

- qubeC3 / qubeC3T Layouteditor wurde um weitere Einstellmöglichkeiten erweitert
- Weitere Einstellmöglichkeiten für Alarmgrenzen wurden geschaffen
- Mediencenter: Erstellen, Ordnen und bessere Übersicht der vorhandenen Medien

## qubePLAN

- Fehler bei der Datumseingabe behoben
- Neues Kursformat NRP verfügbar

## Education Center

- Neben unseren eLearning Lektionen finden Sie hier nun auch die Hyperlinks zur qubeSERIE Anwenderschulung und unserer neuen Wissensdatenbank
- 

## Szenarien

- Mediencenter implementiert. Einfaches Auswählen von vorhandenen Medien aus dem Mediencenter der Organisation.
- Optimierter Freigabe- und Synchronisierungsprozess mit qubeCONTROL führt zu geringeren Ladezeiten

## Kontakt

support@skillqube.com  
+49 6222 – 386265  
knowledge.skillqube.com